

## V、How to use

- 1.Press “ON/OFF” key to turn the massager and red indicator function on or off.
- 2.Press “SPEED+” or “SPEED-” key to increase or decrease the massage intensity level. There are a total of six intensity levels to choose from.
- 3.Press “AUTO” key to use the auto massage mode. There are a total of six massage modes to choose from.
- 4.The massager will automatically turn off after 15 minutes.
- 5.After use, you must turn off the product by pressing the off key.

### Warming

- 1 This product has no timing function, it is recommended not to exceed 15 minutes per use
- 2 When using this product, it may cause interference to the TV, etc. Please keep the distance properly.

## VI、Jugement of Fault

Fault	Reason	Fault Remedy
Motor not working	Plug not fixed up	Fix up the plug well.
	Switch is out of position	Checking the Switch
Working too long	Thermal Protection	Rest for 30 minutes

MADE IN CHINA

soga®

## Infrared Massager

MSAG-L0QQ5B-469147



Please read the handbook carefully before using,especially the part of Caution concerned,and keep the operation book well,please.Thank you.

## Massaging Principle

Massaging is one of outside treatment of Chinese traditional medicine. It can act on the body at specific points by all kinds of technique ,so that it can ajust system physiological and pathology condition. And then you will have a good therapy and strong body.The massaging principle mainly embodies at following points:the changes of system funtion,the ajustment of biological information and unusual correcting the postion of dissection, etc.It is main treatment means or auxiliary measures for some disease,and it is used widely.

According to the theory of Chinese medicine,especially theory of main and collateral channels,it emphasizes that body surface of human is connected with the internal organs by main and collateral channels or points.If the internal organs are something wrong,the disease will reflect on the body surface by main and collateral.And then,we can massage and stimulate the main and collateral or points of body surface,it will make gaseous pulse move well. So,when massaging,the treatment “information” will be sent to the sick internal organs,it will make your internal organs’ function stronger.

From the anatomical knowledge of modern medicine,massaging will promote the natural power of resisting disease,promote absorbing of inflammation exudation,and can dull muscle spasm and ache.

Through practice validation after thousands of years,massaging has a good treatment for inflammation,conglutination,spasm, hyperplasia, degeneration, fibrosis,blood vessels nerve of hurting parenchyma.

## I. Function of Product

- ① To eliminate fatigue and relax muscle because of violent action.
- ② To eliminate ache because of long bad blood circulation.
- ③ To eliminate spasm shoulers because of stiff neck.
- ④ To eliminate ache of tiredness or rheumatism.
- ⑤ Massaging head can massage all the parts of body.

## II、Caution

- ① Do pull out the plug after using or before cleaning it.
- ② Don't use it at moist place.
- ③ When children or the disabled use it,pay more attention to it.
- ④ If the product is something wrong,do not repair it yourself or don't go on using it,you must contact with the seller ASAP.

## III. Technological Parameter

Voltage	220V	
Frequency	50Hz	
Vibration Frequency	High	3200 revolutions/m
	Low	2600 revolutions/m
Power	25 W	

## IV. Chart of body point

There are around 650 points in our body from neck to waist. And every point is connected with some organs or body condition. Massaging can promote blood circulation,promote metabolism and build up the ability agaist disease.Massager will make you have massaging delight and healthy pleasure by easy ways.

